

BioArk® Prestained Protein Marker 180KD

Cat. #: BAPM1125

Product Information

Product Name	Cat.No.	Spec.
Prestained Protein Marker 180KD	BAPM1125-01	1 mL
	BAPM1125-02	5×1 mL
	BAPM1125-05	15×1 mL

Product Description:

This product, Prestained Protein Marker 180KD, is composed of 10 highly pure and pre-stained recombinant proteins. It indicates a molecular weight range of 10-180 kDa in Tris-Glycine gel (~10, ~15, ~25, ~35, ~40, ~55, ~70, ~100, ~130, ~180 kDa), with the 70 kDa markers appearing as pink bands, 10 kDa appearing green band, while the others appear as blue bands. This allows for the dynamic observation of protein electrophoresis or the assessment of protein transfer efficiency. It is suitable as a protein molecular weight standard for SDS-PAGE and Western blot.


Key Features:

(~10, ~15, ~25, ~35, ~40, ~55, ~70, ~100, ~130, ~180 kDa)

Delivery and storage:

Ship with wet ice; Store at -20°C, valid for 12 months.

BAPM1125



Prestained Protein Marker 180KD

For Research Use Only!